MAJOR AF COMP

VENEERING SYSTEM FOR AESTHETIC
PROSTHETIC RESTORATIONS

EN ENGLISH DIRECTION FOR USE
DE DEUTSCH HINWEISE ZUR VERARBEITUNG
FR FRANÇAIS MODE D'EMPLOI
ES ESPAÑOL INFORMACIONES DEL FABRICANTE
IT ITALIANO INFORMAZIONI DEL FABBRICANTE

EN DIRECTION FOR USE

Single-component light-curing material for filling and supplementing of MAJOR AF COMP veneering teeth and of teeth in polymethylmethacrylate material (PMMA) or in composite.

Class IIa MDD - Logo CE 0470

Range of application:

- Fixed restoration:
 - crown, bridges, adhesive bridgework
 - implant-supported super structures, veneers, long term temporary restorations
- Removable restoration:
 - conical and telescopic crowns
 - attachment matrices
 - customizing acrylic teeth in resin for removable restoration

General important points

- MAJOR AF COMP system needs light polymerization with 310-500 nm wave length.
- MAJOR AF COMP Metal Primer must be applied in a thin layer and on sandblasted metal structures.
- Apply MAJOR AF COMP Opäquer - in - a - thin - and homogeneous layer, polymerize, apply a second layer to obtain the correct covering, then polymerize again.
- Each composite layer has to be polymerized separately.
- Avoid composite layers thicker than 1-2 millimeters.
- Do not remove the dispersion layer - should it be removed, it is possible to create a new dispersion layer using the MAJOR AF COMP Bonding.
- Follow the autopolymerization times of MAJOR AF COMP Dual Cem. Working time: 4-5 minutes, the exceeding material must be used as base for further composite (modelling) then polymerize.
- Rotate the plunger clockwise to extrude the paste out from the syringe.
- Do not leave the material exposed to strong light sources. Rotate the plunger anticlockwise to avoid the material leakage. Do not insert the instrument in the syringe!

Direction for use

Preparation of the metallic structure (MAJOR AF COMP Metal Primer).

Once the framework has been casted and finished, treat the structure as usual. Sandblast the surfaces with 110/125 µm, 2/4 bar and use oil free compressed air. After cleaning apply the agent for the metal-resin bond (MAJOR AF COMP Metal Primer). Brush the agent with the minimum thickness coat.

Caution! do not touch the surfaces after cleaning, do not use compressed air or substances which form films.

Apply the primer. The surface should be only slightly moistened with a brush. Allow the surface to dry in the air for 2 minutes. The primer must be completely dry before any opaque is applied.

Applying the opaque (MAJOR AF COMP Opäquer) (MAJOR AF COMP Bondopäquer).

Apply the paste opaque (MAJOR AF COMP Opäquer) in the dental areas; use the desired colour; OL-OM-OD, the three colours can be mixed together to obtain the individual shade. Apply the paste opaque (MAJOR AF COMP Bondopäquer); the OC-PINK colour in the gingival areas. Bondopäquer is a light-curing paste opäquer with an integral adhesive component. Therefore the use of an additional bonding system is not necessary. Apply a thin layer of "Opäquer or Bondopäquer" with a rigid short-haired brush. The complete coverage is obtained applying at least two coats. An non-covering opäquer layer has a negative effect on the chromatic result. The flowability of the opäquer can be optimized by mixing it on a mixing block. Take care that it is applied a particular thin layer in the retention areas. Each layer polymerizes separately. Polymerize for 180 seconds with 310-500 nm wave length. Available shades of the opaque: please see the table in the last page.

Preparation of aesthetic veneering teeth (MAJOR AF COMP Bonding).

Roughen the aesthetic veneering teeth by sandblasting with 110 µm aluminium dioxide at 2/3 bar pressure and clean them with oil free compressed air. Then brush on the sandblasted surface the light-curing bonding for acrylate and composite (MAJOR AF COMP Bonding) and proceed with the prepolymerization to create a new layer of dispersion. Polymerize for 90 seconds with 310-500 nm wave length. A whitish appearance of the pre-polymerized Bonding is the result of an excessive polymerization and it must be removed. Repeat the above mentioned procedure reducing the polymerization time.

Fixing of the aesthetic veneering teeth (MAJOR AF COMP Dual Cem)

MAJOR AF COMP Dual Cem is a dual-component composite with dual polymerization (self hardening and photo-polymerizing) for the fixing of the aesthetic veneering teeth (MAJOR AF COMP) on metal or ceramic structures.

MAJOR AF COMP Dual Cem is applied, according to the selected shade, on the surface of the metal, which is previously opaqued, and on the inner side of the veneering teeth MAJOR AF COMP treated with MAJOR AF COMP Bonding. Veneering teeth are positioned on the metal structure with the help of a silicone key. After the application of the Dual Cem, before the final polymerization it is important to remove the exceeding material and make corrections or fitting. The time limit is of about 4-5 minutes at the room temperature of 23°C. To guarantee a complete polymerization and a good adhesion it is absolutely necessary to harden the product with an irradiation covering an interval of 310-500 nm wave length for at least 180 minutes. Available shades of Dual Cem: please see the table on the last page.

Layering of MAJOR AF COMP/Dentine Flow/Incisal Flow/Gingive Paste/Gingive Flow.

The material is supplied in screw syringes and coded according to its intended use. Available colours: please see the table on the last page. Rotate the plunger of the syringe to extrude the required amount. Apply the material on the surface. Adapt and layer the single material with even thickness not more than 1-2 mm and proceed with polymerization after each layering. Intermediate polymerization for 180 seconds with 310-500 nm wave length. If the individual layer is more than 2 mm thick, prolong the intermediate polymerization time to 8 minutes. The polymerization creates a thin dispersion layer on the surface. This layer is essential for the chemical bonding with the following layers. During the building up this layer must not be removed. Should it be removed, either partially or fully, please refer to the section entitled "Corrections and repairs". After the building up it is advisable to make a cycle of final polymerization. Apply final polymerization for 9 minutes with 310-500 nm wave length.

Finishing

Silicone polishers (wheels, lens-shaped, cylinders), tungsten or diamond burs are suitable for finishing.

Polishing

Polish with goat hair brush, polishing paste and soft cotton discs on hand piece. Finish and polish carefully the surface to obtain a good result and avoid any deposits (nicotine, caffeine, etc.) and colour alterations. **Note:** like the acrylics, the finishing

of composites produces dust. It is advisable to apply the finishing over a dust suction system.

Correction and repair

Roughen the material up to 2 mm above the margin of the area to be repaired or corrected. Apply the MAJOR AF COMP Bonding and pre-polymerize with the light cure unit to create a new dispersion layer. Polymerize for 90 seconds with 310-500 nm wave length. A whitish appearance of the Bonding Fluid is a result of an excessive polymerization and it must be removed. Repeat the above mentioned procedure but reduce the polymerization time. Then proceed with the layering of the composite mass and polymerize.

Possible mistakes and how to avoid them

Opäquer polymerization

- All frame-works have to be degreased and any polishing paste residual should be removed.
- The opäquer must be coated homogeneously. It is better to apply two coats of opäquer, each coat should be separately polymerized.

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- Proper wax-up the framework to ensure the support to the composite materials in any occlusal condition.
- Avoid pre-contacts
- Respect the correct thickness of layers and the polymerization times for the different materials
- Apply retention sphere
- Do not apply too much primer
- Let the primer dry for 2 minutes and then apply immediately the opaque.

Air inclusion

- Rotate the plunger to extrude the paste out from the syringe and take the material with a proper instrument but do not insert it in the syringe.
- Apply sufficient material and spread it.
- Do not mix composite materials but apply them by layers.

Discoloration and plaque

- Check the luminous power of the light curing unit
- Seal the surface by polishing it thoroughly
- Position the frameworks in the light curing unit following the indications
- Apply the composite materials in the correct thickness and follow the polymerization times.

Guarantees

Instructions provided verbally, in written or during practical demonstrations are based on our own experience and should be only considered as guidelines. Our products are subject to continued improvements. We reserve the right to modify them. MAJOR AF COMP complies with EN ISO 10477 standard which does not apply to materials for use in the lateral region which is exposed to masticatory load.

Side effects

Unwanted side-effects of this medical device are extremely rare when the product is regularly used in a proper way. Reaction of the immune system (for instance allergies) or local irritations, however, cannot be completely excluded. Should you learn about unwanted side effects - even if it is doubtful - please contact us.

Contra-indications/reactions

If the patient has hypersensitivities towards a component of this product, we recommend not to use it and to consult the doctor/dentist. Dentists must consider known reactions of the product with other materials already in the patient's mouth before using it.

Note

It this product is used for a special work, kindly supply the dentist with the above information. Please follow the safety data sheets.

Main components of the materials - Inorganic filler (74% by weight or 66% by volume) - Glass filler (medium granulometry: 0.7 µm) - Pyrogenic silica (medium granulometry: 0.04 µm) Monomers (25% by weight) - Urethane dimethacrylate - Butanediol dimethacrylate - Bis GMA - Additives (1% by weight) - Catalysts, stabilizers, pigments

Working times of materials:

1-3 minutes, depending on the light conditions

Depth of polymerization

The depth of polymerization is 2.0 mm. To ensure optimum physical characteristics, no layer should be more than 2.0 mm thick. Use light curing units with emission spectrum of 310-500 nm. The requested physical characteristics can only be obtained with perfect lamps. It is necessary to regularly check the light intensity following the supplier's information.

Storage

To guarantee the correct shelf life do not store this material above 25°C. If MAJOR AF COMP is kept in the refrigerator (in case of high outside temperature), it must be removed in good time before use to resume its modelling properties. Close the syringe immediately after use. Protect the syringe against direct sun light. Turn the plunger counter-clockwise to avoid the material to leak out.

Shelf-life

The shelf life is printed on the label of each syringe. Do not use the product after the expiry date.

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