VILLACRYL T MONO/MULTI INSTRUCTIONS FOR USE - EN

Villacryl DISC T Mono/ Multi – PMMA disc for milling dental arches and teeth in removable dentures

Product description

Villacryl DISC T (Mono/ Multi) are PMMA discs designed for CAD/CAM technology to produce dental arches in removable dentures.

Indications for use: Complete dentures (edentulism), partial removable dentures.

Clinical benefits: Restoration of chewing and biting functions lost due to partial or total loss of natural teeth.

Target group(s) of patients and intended users of the products:

Users: Professionals in the dental technology industry - certified dental technicians making dental prostheses using CAD-CAM technology.

Patients: people of all ages and genders, with partial or complete loss of teeth, eligible for prosthetic treatment whereby removable dentures, including children, adolescents, and adults.

Contraindications

Do not use the product for people with a diagnosed allergy to the ingredients of the product. In such situations, use a different material, suitable for the patient.

Adverse effects

Dental technicians

Irritation, redness or hypersensitivity symptoms including anaphylactic shock if allergic to any of the ingredients of the product.

Patients

Irritation, inflammation and allergic reactions of the oral cavity caused by fungi and bacteria adhering to the denture.

Product form

Disc 16 or 20 mm thick and 98.5 mm in diameter, in various colour hues similar to natural teeth.

Forms of delivery: PMMA disc wrapped in a cardboard box and protective foam.

Product composition: methyl methacrylate polymer, pigments.

Compatible products

For use with acrylic materials for denture plates, cured or milled using CAD CAM technology, and materials for bonding artificial teeth into a milled denture plate.

Safety instructions

Warnings

When coming into contact with the technician's body, dust generated during milling can irritate the skin, cause an allergic skin reaction or irritate the respiratory tract.

If dust particles enter the eyes, mechanical irritation of the eyes may result

Precautions

Avoid inhaling dust generated during milling.

Wear protective gloves/protective clothing/eye protection/face protection.

Do not eat or smoke while working with the material.

Use ventilation in the workplace.

If dust comes into contact with eyes, flush using water and consult a doctor.

Wash areas of skin that have come into contact with the product ingredient thoroughly using soap and water.

Information to be provided to patients / denture users

If symptoms of allergic reactions or irritation occur, discontinue use of the denture immediately and consult a doctor.

The strong oxidising agents in cleaning products and mouthwashes can change the colour of the teeth during use.

After each meal, the denture should be washed with a soft toothbrush under cool running water. Unscented soap, various denture cleaners and denture care products that do not contain abrasive substances can be used for cleaning.

Teeth made from the material are not resistant to prolonged exposure to alcoholic beverages.

When hygiene of the denture is inadequate, substances contained in the liquids consumed (coffee, tea, coca-cola, etc.) can cause deposits and discolouration on the surface of the teeth.

Aggressive chemicals used in denture cleaners can change the colour of the teeth in the denture

or cause them to lose their original physical properties.

Any damage to the teeth in the denture should be reported to the dentist. Repairs can only be carried out by a dental technician.

If a piece of damaged tooth is swallowed, consult a doctor. The material used to create teeth does not show in X-ray examination.

Storage conditions

Store the product in the original packaging in a dark room, at a temperature of 5-25 °C, away from unauthorized persons.

Waste disposal

Dispose of in accordance with current business waste management regulations.

Cardboard packaging, labels, instructions: recyclable, dispose of as non-hazardous separated municipal waste - paper.

Acrylic disc: non-hazardous segregated municipal waste - plastic.

Directions for use (step-by-step instructions)

Preliminary activities:

Read the instructions carefully before using the product

Before making the prosthesis, ensure that the patient for whom the prosthesis is being made is not allergic to any of the material's ingredients.

NOTE: In order to choose the right tooth colour, it is recommended to use a colour matching chart (colour chart). Colour should be chosen individually and if in doubt, do your own test. Before using the disc, read the disc label carefully to select the correct disc height and hue. To design a dental arch or a single tooth, use appropriate CAD software. When designing, the minimum wall (tooth) thickness must be 0.6 mm. Failure to observe the minimum thickness may result in damage (e.g. cracking). Villacryl DISC T (Mono/ Multi) can be milled in CAD/CAM milling systems using standard 98.5 mm chuck. Only the milling tools recommended for this system should be used to process the material. Note that the material may overheat during dry milling. This can damage the material. In such case, we recommend wet milling.

Bonding teeth to the denture plate

Artificial teeth made from Villacryl DISC T (Mono/ Multi) and the denture plate should be bonded using PMMA self-curing resin (Villacryl S) or photocurable materials, according to the manufacturer's instructions. The mucosal surface in artificial teeth should be roughened and rubbed twice with acrylic monomer. In addition, bond strength can be increased by applying Villacryl Bond material to the tooth surface and milled denture plate.

Finishing and polishing of the prosthesis

After milling, separate the milled teeth from the disc abutments using a suitable cutter. To finish or adjust the shape, use a cross-cut cutter. Particular attention should be paid to the interdental spaces.

Polishing

Apply standard polishing process, using polishing tools suitable for dental acrylic materials.

Adjustments and repairs of damaged dentures

Teeth made from Villacryl DISC T (Mono/ Multi) can be repaired or adjusted using conventional self-polymerising acrylic resin, after you first wet the item to be adjusted or repaired using liquid.

Packaging and transportation of the denture to the dental surgery

Finished prosthesis should be protected from mechanical damage and high temperatures.

Hygiene, sterilisation, disinfection

Product

The product is designed for the fabrication of a single prosthesis. It is not supplied sterile, nor is it required to be sterilised or disinfected before use. Apply the principles of good laboratory practice in dental laboratories to protect the product from contamination.

Prosthesis made using the product

There is a risk of cross-infection. The prosthesis should be washed and disinfected by the dentist before being placed in the patient's mouth for the first time.

For disinfecting dentures made of the product, it is recommended to use disinfectants intended by their manufacturers for disinfecting dental prostheses and to follow the manufacturer's recommended procedures.

Agents based on at least 65% ethyl or propyl alcohol can be used. The denture should be sprayed with a disinfectant solution and then gently wiped dry with a paper tissue. Ready-to-use disposable wipes soaked in a disinfectant solution can also be used. In addition, the denture should be washed with

warm water and soap. Communicate the relevant information to the dentist.

Additional information

Reporting serious incidents: Any serious incidents involving the products must be reported to the manufacturer and to the competent authority in accordance with local regulations.

Data necessary for the identification of the product and package contents

Reference number, product name, colour, disc thickness, batch number and expiry date - printed on the product packaging.

The expiry date refers to the unopened packaging properly stored. Do not use the product after the expiry date.

Packaging contents

VILLACRYL DISC T (Mono/ Multi)

Colour	Available thickness in mm
A1	16, 20
A2	16, 20
A3	16, 20
B1	16, 20

Prevention of problems

See website www.everall7.pl.

Please send any comments on our products to our address.

For details on all Everall7 materials, please call: +48 602 701 001

Information provided in any way, including demonstrations, will not invalidate the instructions for use. Users are required to make sure that the product is suitable for the intended use.

The manufacturer shall not be held liable for damage, including to third parties, resulting from a failure to follow the instructions or unsuitability for use. The manufacturer's liability is in any case limited to the value of the product supplied.

MD

The symbol means a medical device as defined in EU Regulation 2017/745.

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