Gebrauchsanweisung Dentify XACT

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Instruction for use Identify XACT

Identify XACT is a phosphate-bonded precision investment especially developed for the production of crowns and bridges. It is suitable for all Dental alloys (precious and non precious). By respecting the following working details, time and temperature specifications, you will achieve precise results.

13. Application / Advantages

Precision investment for crown and bridge as well as for implant technique. Can be used in conventional or Speed casting method.

- Easy to control the Expansion
- Fine and dense consistency
- Smooth surfaces and precise fitting

Identify XACT

Art-No.: 731-0006

Content: 38 Bags of 160g

Identify ONE Universal liquid

Art-No.: 710-0001

Content: 1 Litre Liquid

14. Working place

The casting results depend on the cleanness of the user. Even small residues on the working units can affect the casting results. Please always use special Spatula and Mixing Bowl when mixing the investment and separate instruments when working with stone and plaster.

15. Storage

The powder and the liquid of the investment must be stored at a room temperature (20-23°C). The Identify ONE Universal liquid must be protected against freeze and should not be placed or stored in a place where the temperature is below than 5°C, that may cause damage of the product and will not be suitable for use anymore.

16. Danger / Health warning

Investment contains Quartz! Please pay attention not to breath dust during work. It might cause Lungs diseases (for ex. Lung cancer, Silicosis). Use a suitable dust mask and avoid any dust inhalation when emptying the investment powder in the mixing bowl. Remove dust residues from the working place with a wet cloth. Wash the empty bags before throwing them away. Use a filter and wet the investment ring before deflasking to avoid always the dust inhalation.

17. Technical Data

Working time: 6 – 7 minutes
Working temperature: 20 – 23°C Powder and Liquid
Mixing ration: 160g Powder : 42 ml Liquid

18. Mixing ratio

Identify XACT should be mixed with Identify ONE Universal liquid

The below recommendations are guide lines based on our own test results. Too many factors, such as humidity, room temperature etc. can influence the results.

Metal (160g : 42ml)

	Crowns & Bridges of Wax		Post & Core Inlays	
	Liquid	dist. Water	Liquid	dist. Water
Non precious alloy	38 ml / 90%	4 ml / 10%		
Precious alloy with higher Gold content (70-80% Au)	20 ml / 47%	22 ml / 53%	18 ml / 42%	24 ml / 58%
Precious alloy with reduced Gold content (55-65% Au)	21 ml / 50%	21 ml / 50%	18,5 ml / 44%	23,5 ml / 56%
Palladium alloy	28,5 ml / 68%	13,5 ml / 32%		

Hint:

By increasing or decreasing the amount of the distillate water to the Identify ONE Universal liquid, the expansion can be controlled. The total amount of the liquid quantity 160g : 42 ml shall not be modified.

In general, more you add Identify ONE Universal liquid, more will be the expansion. More you add distillate water, less will be the expansion. The expansion can be slightly influenced by modifying the total amount of the liquid (max. 4ml). In this case the following rule will be taken into consideration, a thicker mix will lead to a bigger expansion and a thinner mix will lead to a lower expansion.

In case of very special work case, such as thin coping implant constructions, and in order to reach the desired stability, we recommend you to reduce the total amount of the mixed-up liquid. But also in order to preserve on the exact fitting, the amount of the ONE universal liquid should be also reduced in such cases.

19. Preparation

After positioning the casting sprues and made all the necessary steps to prepare the model for imbedding, you can apply a thin layer of debubblizer on the wax modellation to reduce the tension and avoid the creation of air bubbles.

20. Processing / Mixing

The working process should be done at room temperature 20- 23°C (Warm shorten the working time). Take into consideration that the mixing bowls for the corresponding powder are clean of any residues. They should not be completely dry and not to be used for stone or investment containing plasters, as these can influence the quality of the results.

Weigh the powder and fill it in the mixing bowl. Add the mixed up total amount of liquid. Start measuring the time and count exactly 20 minutes then the mould is ready. Mix well first by hand for 15 seconds till the consistency becomes even, then under vacuum for 45 seconds. (Tip: The mixing unit can also influence the quality of your results, that's why we recommend you to control it from time to time. A bad mixing unit leads to fitting tolerance and air bubbles in the casting objects. Mixing speed should be adjusted at 250 r/min. After leaving the investment for 10 further seconds in the vacuum unit you can start filling in the investment under low vibration. At a temperature of 23°C the working time could be 6 min. since the beginning of the mixing process.

Special tips for resin modellation:

Leave the resin modellation for further 10 minutes in a pressure pot at 45 – 55°C for polymerisation. Place the muffle after 20 minutes in a preheated furnace at 360°C to burn out the resin. Holding time is 30 Min. Heat up to the desired final temperature without further holding time or rather put it in another preheated furnace already heated at final temperature.

21. Heating

Place the mould into the preheating furnace with the cone face down onto a punched or coarse milled plate.

a. Conventional heating process:

Place the mould into a cold furnace, soonest after 20 minutes since the beginning of the mixing process. Start heating up with a heating rate of 20°C/Min. Holding time is not necessary only when working with resin. Hold at final temperature for at least 45 minutes according to the alloy instruction for use.

b. Speed heating process

Place the mould after 20 minutes, since the beginning of the mixing process, into the hot furnace for 45 minutes at a maximum temperature of 850°C. Heat up further to the end temperature if the alloy used requires that, please refer to the alloy instruction for use.

The holding time depends on the size and number of the muffles. In case of several, then extend the holding time 10 minutes for each.

22. Casting

After a holding time of at least 45 minutes at final temperature the casting can be started according to the instruction for use of the alloy. The casting should be achieved immediately after taking out the muffle from the heating furnace.

23. Cooling

Cool down the muffle slowly under room temperature, placing the casting funnel upwards.

24. Warranty

We guarantee the quality of our product. The given technical instructions are based on our test results. These can be only guaranteed if the processing has been carried out as mentioned above. Every user is responsible by himself for the correct appliance of these rules and given data into the practical field. We, are not liable for any faulty results as Dently GmbH has no influence on the working procedures done by the user. The given data corresponds to the actual technical stand. Our products are permanently subjected to further development. For this reason we keep the right for any modifications in the production or components. Claims on any damages are only related to the quality of the products.