## PL- Zapobieganie problemom

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## EN- Problem prevention

Problem	Potential impact on patient / device	Presumed cause	Recommended action
<b>An allergic reaction of the patient</b>	-	<i>The patient is sensitive to allergenic compounds contained in the product</i>	<i>Use other material suitable for the patient</i>
	-	<i>Increased level of residual monomer</i>	<i>Strictly follow the powder/liquid ratio</i>  <i>Use recommended way of reduction residual monomer level in the IFU.</i>
<b>Porous structure in the finished prosthesis</b>	<i>Increased residual monomer level - may cause irritation and allergic reaction</i>  <i>Reduction of the mechanical properties of the prosthesis, breaking of the prosthesis during use, veneered surface peeling off the metal structure</i>	<i>Incorrect powder/liquid ratio</i>	<i>Strictly follow the powder/liquid ratio</i>
		<i>Incorrect process of pressing the acrylic dough in the flask under the hydraulic press</i>	<i>Use the pressure and pressing time specified in the IFU.</i>
		<i>Incorrect polymerization process</i>	<i>Use the recommended process parameters.</i>
		<i>Not enough of acrylic material in the flask</i>	<i>10-15% more acrylic mass should be put in the flask. The excess material coming out of the can must be visible.</i>
<b>Reduction of the mechanical properties</b>	<i>Breaking out veneered surfaces during usage</i>	<i>Incorrect powder/liquid ratio</i>	<i>Strictly follow the powder/liquid ratio</i>
		<i>Incorrect polymerization process</i>	<i>Use the recommended process parameters.</i>
<b>Rough surface of the veneered surface</b>	<i>Increased affinity for the accumulation of fungi and microbial flora - possible irritation and microbial infections</i>	<i>Incorrect isolation of stone model</i>	<i>Check the type of isolator used. A film-forming insulator, eg Izo-sol, should be used.</i>
<b>Incorrect color of veneered surface</b>	<i>The esthetics of the prosthesis is not acceptable</i>	<i>Incorrect powder/liquid ratio</i>	<i>Strictly follow the powder/liquid ratio</i>
		<i>Badly shaped tooth. Too much or too little enamel layer compared to the dentine.</i>	<i>Shape and fabricate a tooth from scratch.</i>
		<i>Badly matched enamel and dentin colors</i>	<i>Use the recommended enamel and dentin shades listed in the IFU table</i>
<b>Change in color and discoloration of the denture during use</b>	<i>The esthetics of the prosthesis is not acceptable</i>	<i>Soaking the technician's hands with water / using moisturizing creams when processing the material</i>	<i>Use polyethylene HDPE gloves</i>
		<i>Sludge from smoking, drinking coffee, dyes contained in food or spices used in food, medications</i>	<i>Maintain proper oral hygiene.</i>
<b>The metal structure is visible under the veneered tooth.</b>	<i>The esthetics of the crown is not acceptable</i>	<i>Too thin acrylic layer</i>	<i>Shape and fabricate a tooth from scratch.</i>
		<i>No opaque layer on the metal structure to be veneered</i>	<i>Apply Villacryl Opaker in the appropriate shade to the veneered tooth.</i>

<b>Breaking out veneered surface on metal structure</b>	<i>Crown damage. Loss of intended properties and functions.</i>	<i>Poorly prepared metal surface for veneering.</i>	Make mechanical retention on the metal surface for better bonding with the material.
		<i>Incorrect polymerization process</i>	<i>Use the recommended process parameters.</i>
<b>The gap between the veneered and metal structure is too large</b>	<i>Increased affinity for the accumulation of fungi and microbial flora - possible irritation and microbial infections</i>	<i>Polymerization shrinkage</i>	<i>Strictly follow the powder/liquid ratio</i>  <i>Use the recommended process parameters.</i>